


The Bronchus Through the Bronchovideoscope

Normal and Abnormal Findings



A 59-year-old man. Findings of the right upper lobe segmental bifurcation. B³ is obstructed by a lid-like glossy and flat tumor. There are no findings of invasion to surrounding mucosa. The appearance suggests a benign tumor or a low malignancy epithelial tumor. It was surgically removed by a high-frequency snare and was diagnosed as bronchial lipoma.

An 81-year-old woman. A nodular tumor covered by necrosis obstructs the left basal bronchus. It is the proximal tip of a tumor originating in the left lower lobe, proliferating endobronchially as a polyp. The upper left of the photograph shows a part of the left B⁶. Poorly differentiated adenocarcinoma.

Squamous cell carcinoma originating in the right lower lobe bronchus in a 66-year-old man. A nodular protrusion type tumor projects from the membranous portion of the right main bronchus. The surface of the tumor bleeds easily and irregular depressions and protrusions are recognized. The membranous portion which continues to the base of the tumor, contains invasion of the extra muscular layer, and the mucosal longitudinal folds are elevated.

A 56-year-old man. Left B³_{DEC} is obstructed by a nodularly protruding tumor. The longitudinal folds that continue to the tumor are compressed and thickened, indicating intramural invasion. This suggests that the tumor originated peripherally and invaded from beyond the bronchial wall, with endobronchial growth in part. Squamous cell carcinoma.

A 68-year-old man was referred for common cold-like symptoms and middle lobe syndrome was diagnosed clinically. This shows the carina viewed from the trachea. The lesion extended widely as multiple nodular protrusions from the anterior wall of the trachea to the right lower lobe bronchus. Biopsy yielded a diagnosis of primary tracheo-bronchial amyloidosis.


Disclaimer: Any content or information ("Content") presented herein is illustrative in nature and does not guarantee or represent specific information, outcomes, or results. Olympus Corporation, its subsidiaries, affiliates, directors, officers, employees, agents, and representatives (collectively "Olympus") does not represent to or warrant the accuracy or applicability of the Content. Under no circumstances shall Olympus be liable for any costs, expenses, losses, claims, liabilities, or other damages (whether direct, indirect, special, incidental, consequential, or otherwise) that may arise from, or be incurred in connection with, the Content or any use thereof.

OLYMPUS®

Your Vision. Our Future


Kenkichi Oho, M.D., Ph. D.
Kazuyo Matsukawa, M.D., Ph. D.
Shinkawabashi Hospital, Kawasaki, Japan

Right lateral



Upper Lobe	B¹: R. apicalis	<ul style="list-style-type: none"> a. Rm. apicalis proprius b. Rm. subapicalis ventralis
	B²: R. lobi superioris dorsalis	<ul style="list-style-type: none"> a. Rm. subapicalis dorsalis b. Rm. lobi sup. horizontalis
	B³: R. lobi superioris ventralis	<ul style="list-style-type: none"> a. Rm. lobi sup. ventr. lateralis b. Rm. lobi sup. ventr. medialis
Middle Lobe	B⁴: R. medius lateralis	<ul style="list-style-type: none"> a. Rm. lateralis b. Rm. medialis
	B⁵: R. medius medialis	<ul style="list-style-type: none"> a. Rm. superior b. Rm. inferior
Lower Lobe	B⁶: R. lobi inferioris superior	<ul style="list-style-type: none"> a. Rm. superior b. Rm. lateralis c. Rm. medialis
	B*: R. lobi inferioris subsuperior	
	B⁷: R. mediobasalis	<ul style="list-style-type: none"> a. Rm. ventralis b. Rm. dorsalis
	B⁸: R. ventrobasalis	<ul style="list-style-type: none"> a. Rm. lateralis b. Rm. basalis
	B⁹: R. laterobasalis	<ul style="list-style-type: none"> a. Rm. lateralis b. Rm. basalis
	B¹⁰: R. dorsobasalis	<ul style="list-style-type: none"> a. Rm. dorsalis b. Rm. lateralis c. Rm. medialis

Left lateral



Upper Lobe	(Upper division bronchus)	B¹+²: R. apico-dorsalis	<ul style="list-style-type: none"> a. Rm. apicalis b. Rm. subapicalis dorsalis c. Rm. lobi sup. horizontalis
		B³: R. lobi superioris ventralis	<ul style="list-style-type: none"> a. Rm. lobi sup. ventr. lateralis b. Rm. lobi sup. ventr. medialis c. Rm. lobi sup. ventr. superior
	R. lingualis	B⁴: R. lingualis superior	<ul style="list-style-type: none"> a. Rm. lateralis b. Rm. ventralis
Lower Lobe		B⁵: R. lingualis inferior	<ul style="list-style-type: none"> a. Rm. superior b. Rm. inferior
		B⁶: R. lobi inferioris superior	<ul style="list-style-type: none"> a. Rm. superior b. Rm. lateralis c. Rm. medialis
		B*: R. lobi inferioris subsuperior	
		B⁸: R. ventrobasalis	<ul style="list-style-type: none"> a. Rm. lateralis b. Rm. basalis
		B⁹: R. laterobasalis	<ul style="list-style-type: none"> a. Rm. lateralis b. Rm. basalis
		B¹⁰: R. dorsobasalis	<ul style="list-style-type: none"> a. Rm. dorsalis b. Rm. lateralis c. Rm. medialis