

EndoClot® Polysaccharide Hemostatic System (PHS)

Quick Reference Guide

CAUTION: This guide is only a summary of the instructions for use for the EndoClot Hemostatic System. Be sure to follow the instruction manual included in the device packaging.

Getting Started



1. Ensure that the EndoClot Air Compressor is powered on, the green light is on "H" and that air is exiting the white tube.

2. Open EndoClot Applicator. Peel back cover on package.



3. Open EndoClot Absorbed Modified Polyper (AMP®) particles. Ensure that the correct volume (3g or 5g) AMP particles are chosen.

4. Connect EndoClot Applicator to Air Compressor Tubing. Connect the blue gas filter to the Air Compressor white tube.

5. Attach the AMP particle dispenser firmly onto the end of the gas/powder mixing chamber.

EndoClot® PHS

Quick Reference Guide

CAUTION: This guide is only a summary of the instructions for use for the EndoClot Hemostatic System. Be sure to follow the instruction manual included in the device packaging.

Insert catheter into endoscope



6. Insert EndoClot applicator catheter into endoscope channel until catheter exits tip of endoscope. Avoid bending catheter while inserting.
Position tip of catheter 1-2 cm from the bleeding site.

CAUTION Direct contact of the catheter tip with any fluid must be avoided to prevent blockage. If the catheter is occluded, replace it with a new one, and then repeat the process.



Begin flow of particle



NOTE: Angles greater than 60 degrees during this step may cause occlusion of the product in the particle dispenser.



NOTE: If bleeding continues, irrigate excess AMP particles with a saline solution, then re-apply.

7. Gently tilt the particle dispenser (at a 45-degree to 60-degree angle) to begin the flow of powder into the mixing chamber. This is the "on" position for delivering EndoClot PHS to the bleeding site.

Ensure the tubing is facing away from you.

8. Deliver powder to bleeding site.

While particle dispenser is in the "on" position (upside down with air tubing facing away), begin tapping gently on the back of the particle dispenser. Spray a liberal amount of AMP® particles on the site.

The AMP particles can be irrigated at any time to confirm hemostasis.

9. When hemostasis is achieved, return the particle dispenser to an upright position to stop the flow of particles.

Once hemostasis is achieved, remove excess AMP particles with irrigation.

Additional treatment modalities may be applied at this point.

EndoClot® PHS

Quick Reference Guide

CAUTION: This guide is only a summary of the instructions for use for the EndoClot PHS. Be sure to follow the instruction manual included in the device packaging.

Insert catheter into endoscope

Application Techniques:


- Removing excess blood from the site of bleeding prior to application is essential to achieve optimal hemostatic efficacy.
- Hold the AMP particle dispenser in an upright position to avoid spillage when inserting the catheter into the endoscope.
- To avoid occlusion, keep the catheter tip 1-2 cm away from the lesion at all times.
- Direct contact of the catheter tip with any fluid must be avoided to prevent blockage caused by reverse flow. If the catheter is occluded, replace it with a new one, and then repeat the process.
- Follow your standard guidelines for reprocessing of the endoscope.

Troubleshooting

Air Compressor won't turn on

Ensure that power is connected to the air compressor and that the switch on bottom of unit is turned to AC Power

If battery power is utilized, ensure that the batteries are fully charged and the switch on bottom of unit is turned to Battery Power



After using pump with battery power, switch unit back to AC power to preserve battery life

Olympus is a registered trademark of Olympus Corporation, Olympus America Inc., and/or their affiliates. | Manufactured by EndoClot Plus Co., Ltd. Medical devices listed may not be available for sale in all countries.

OLYMPUS

OLYMPUS AMERICA INC.
3500 Corporate Parkway, PO Box 610, Center Valley, PA 18034

For more information, contact your Olympus sales representative, or call 800-848-9024.
www.medical.olympusamerica.com
©2022 Olympus America Inc. All rights reserved.
Printed in USA OAIET0322QRG45234