Eliminate Manual Cleaning Steps with the OER-Elite™ Endoscope Reprocessor

Endoscope Cleaning Guide for EVIS EXERA II (180) & EVIS EXERA III (190) GI Endoscopes

WARNING

When reprocessing the TJF-Q180V and TJF-Q190V duodenoscopes using the OER-Elite endoscope reprocessor, conduct precleaning and manual cleaning as detailed in each duodenoscope's reprocessing manual. This guide is only a summary of the steps necessary to properly reprocess your endoscope. Be sure to follow the detailed steps outlined in the ENDOSCOPE REPROCESSING MANUAL that was included with your endoscope when purchased.

PRECLEANING – Immediately After Use

1

Wipe down insertion tube with a water-soaked cloth.

2

Aspirate water, then air through the suction channel.

3

Attach AW channel-cleaning adapter. Flush water, then air, through air/water channels.

4

Flush water, then air, into elevator-wire channel/auxiliary-water channel. (Alternatively, the auxiliary-water channel may be flushed using the OFP pump.)

For endoscopes with elevators or auxiliary-water feeding only:

5

Disconnect all detachable parts.

LEAKAGE TESTING

6

Attach water-resistant cap. This step is omitted for 190 endoscopes. Transport to the reprocessing room in a covered container. **NOTE:** Ensure that the ETO cap is removed prior to leak testing 190 generation scopes.

7

CAUTION: If a leak is detected, the endoscope must be repaired. Follow manufacturer's instructions.

Ensure the ETO valve and the MB-155 are free of water drops. Detach auxiliary water tube and attach leakage tester. Turn on air source and confirm expansion of the bending section. Immerse entire endoscope in clean water. Perform leakage test. Angulate tip during test. Remove from water. Turn off air source and detach leakage tester from air source. After the tip has deflated, detach leakage tester from endoscope.

MANUAL CLEANING

8

While immersed:

Immerse in freshly prepared detergent solution. Clean all external surfaces. Brush endoscope distal tip (and forceps-elevator, if applicable).

9

While immersed:

Brush the insertion tube portion of the suction channel. Repeat until all debris is removed.

10

While immersed:

Brush the universal cord portion of the suction channel. Repeat until all debris is removed.

11

While immersed:

Brush the suction cylinder and instrument channel port. Repeat until all debris is removed.

12

Suction Cleaning Adapter
Distal tip

Aspirate detergent solution into suction channel.

13

Injection Tube
Channel Plug

Attach injection tube and channel plug. Flush detergent solution into the air and water channels. Flush the forceps-elevator, if applicable.

14

While immersed:
For endoscopes with elevators or auxiliary-water feeding only:

Washing Tube

Flush detergent solution into auxiliary-water / elevator-wire channel.

15

Wipe all external surfaces of endscope and cleaning accessories, then soak in detergent solution for the time recommended by the detergent manufacturer.

Can be omitted if reprocessing with the OER-Elite endoscope reprocessor

RINSING

16

Immerse the entire instrument in clean water and gently agitate to rinse.

17

While immersed:

Auxiliary-
Injection Tube
Channel Plug

Flush clean water through all channels (including auxiliary-water / elevator-wire, if applicable). Remove from water, then flush air through all channels.

18

Use a soft, lint-free cloth to remove excess moisture from endoscope and cleaning accessories in preparation for disinfection.

IMPORTANT

- Wear all appropriate personal protective equipment.
- Meticulous cleaning is essential for effective disinfection/sterilization.
- Be sure to reprocess all removable parts (e.g., valves) according to the reprocessing manual.
- Visually inspect the equipment after cleaning. If debris remains, repeat the procedure.
- Inspect the MB-155 (Leakage Test Connector) for signs of damage prior to use.

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