

Using DualKnife J's jet function

Movie 4 : Tips for mucosa lifting

In cases with fibrosis, the knife's jet function is particularly useful. Fig. 10 shows a mid-procedure phase in which the fibrotic region has been exposed by digging on both sides. Because secure preservation of the dissection space is desirable at this stage, the usual approach is to switch devices and insert a local injection needle. However, changing devices takes time, and it is not uncommon for the field of view to change during that interval. Instead, using the DualKnife J's jet function allows you to proceed rapidly with dissection while delivering fluid. When local injection is performed with a needle, the submucosal cushion may dissipate when the device is switched to the knife. The jet function, on the other hand, makes it possible to start dissection at the exact moment the submucosal layer is elevated.


The procedure to use the jet function is simple: 1) Retract the tip; 2) Lightly touch the target area with the sheath tip to deliver the liquid press the coagulation wave pedal on the electrosurgical unit; 3) Project the knife tip; 4) Start the jet function through the knife.

This enables effective elevation of the submucosal layer. Efficient delivery is also possible by connecting a water-jet pump to the knife. If delivering saline, you can attach a 5–10 cc syringe externally to the knife and push by hand (Fig. 11).


In this case, it would also be acceptable to use a traction device at the initial stage. We will show an example using traction in another case: an LST that appears to enter the ileocecal valve. Although the terminal ileum side of the lesion is difficult to see, applying traction after incision makes it easier to identify the resection line. Applying a second traction device in a different direction can further facilitate resection (Fig. 12).

Be careful when using traction. Traction in the wrong direction will not yield the expected result, and excessive traction can tear the mucosal flap. In addition, the muscular layer is often pulled by traction, which can lead to perforation if the dissection line is mistaken. Therefore, it is necessary to pay close attention to the resection line (Fig. 13).

Fig. 10



a: The arrowhead indicates the dissection space. Because fibrosis is present immediately to the left, it is desirable to expand the space before proceeding with dissection.




b: Jet function from the tip of DualKnife J expands the dissection space indicated by the arrowhead.


Fig. 11




Fig. 12



a: LST partially entering the ileocecal valve.




b: Incising the terminal ileus side and then lifting with a traction device allows visual recognition of the resection line.



c: Using another traction device facilitates resection further.

Fig. 13



Because traction also pulls the muscular layer, perforation would occur if dissection were performed along the blue line.

Reference materials: 1) Miura Y, et al. Gastrointest Endosc. 2016;83(2):457-8. 2) Yahagi N, et al. Endoscopy. 2017;49(10):E227-E228.


Any content or information ("Content") presented herein is illustrative in nature and does not guarantee or represent specific information, outcomes, or results. Olympus Medical Systems Corp. and its parents, subsidiaries, affiliates, directors, officers, employees, agents, and representatives (collectively "Olympus") does not represent to or warrant the accuracy or applicability of the Content. Under no circumstances shall Olympus be liable for any costs, expenses, losses, claims, liabilities, or other damages (whether direct, indirect, special, incidental, consequential, or otherwise) that may arise from, or be incurred in connection with, the Content or any use thereof.

**OLYMPUS MEDICAL SYSTEMS CORP.**  
Olympus Medical Systems Corp. 2951 Ishikawa-machi, Hachioji-shi, Tokyo 192-8507, Japan

Trademarks, product names, logos, or trade names used in this document are generally registered trademarks or trademarks of each company.

The clinical movies introduced in this report can be found at this link.





# Colorectal ESD using DualKnife J

Project associate professor, Next-Generation Endoscopic Computer Vision/Gastroenterology, University of Tokyo Hospital

## A/Prof. Yosuke TSUJI

Biography

2004 Mar. Graduated from Faculty of Medicine, The University of Tokyo  
Apr. Resident at JR Tokyo General Hospital  
2006 Apr. Dept. of Gastroenterology, JR Tokyo General Hospital  
2007 Apr. Dept. of Gastroenterology, The University of Tokyo Hospital  
2008 Apr. Dept. of Gastroenterology, NTT Medical Center Tokyo  
2013 Apr. Dept. of Gastroenterology, The University of Tokyo Hospital  
2014 Mar. Master's degree at Graduate School of Medicine, University of Tokyo

2014 Apr. Dept. of Endoscopy and Endoscopic Surgery, University of Tokyo Hospital  
2015 Feb. Assistant Prof. at Dept. of Gastroenterology, University of Tokyo Hospital  
2016 Jan. Part-time teacher at Dept. of Internal Medicine III (Gastroenterology), Kyorin University  
2022 Apr. Project Lecturer, Dept. of Gastroenterology, University of Tokyo Hospital  
2023 Jan. Project Associate Professor, Next-Generation Endoscopic Computer Vision/Gastroenterology, University of Tokyo Hospital

Introduction to the Hospital The University of Tokyo Hospital

Location: Hongo 7-3-1, Bunkyo-ku, Tokyo, Japan / Number of beds: 1,157 (as of FY2024)  
The Department of Gastroenterology of the hospital is composed of three groups: GI tract diseases, liver diseases and biliary & pancreatic diseases. All three groups are actively involved in both clinical and research domains. The department also makes active efforts in endoscopic treatment in close collaboration with other departments – mainly surgical – that deal with high-difficulty, high-risk cases, most of which are referred from other hospitals. At present, the department performs ESD in about 100 esophageal cases, about 160 gastric cases, and about 200 colorectal cases every year. In recent years, it has also focused on endoscopic full-thickness resection (EFTR) of gastric submucosal tumors, endoscopic duodenal treatment, and hypopharyngeal ESD jointly performed with the ENT Dept.



Endoscopic submucosal dissection (ESD) was initially developed for gastric lesions in the late 1990s and was subsequently adapted for colorectal lesions. More than a quarter of a century has passed since the technique was first developed. Over this time, exciting new devices have been developed and methodologies have also evolved. Ideally, each physician should use the procedures and devices that suit them best. Needle-type knives are used universally in current colorectal ESD in Japan. This article discusses strategies for colorectal ESD in actual cases using needle-type knives.

Knives for use in colorectal ESD

Before presenting cases, I will discuss the incision device I use. As for "Which knife is best?", the optimal answer is to select the one that best aligns with your preferences, based on an understanding of the characteristics of each knife. Because the muscularis propria is thin, colorectal ESD requires precise incision and dissection. Needle-type knives are very useful from this standpoint.

The DualKnife J has a knob-shaped tip. With a tip diameter of 0.65 mm and a knife diameter of 0.4 mm, it can hook incised tissue without slipping and allows steady insertion of the knife tip even into narrow gaps, thanks to its moderate diameter. In addition, it has a jet function which allows dissection while delivering fluid to the submucosal layer.



Product name	DualKnife J	
Model	KD-655L	KD-655U
Working length	1,650 mm	2,300 mm
Cutting knife length when extended	2 mm	1.5 mm
Needle tip length when retracted	0.1 mm	
Knife diameter	0.4 mm	
Needle tip diameter	0.65 mm	



## [Case 1] Rectal Ra LST-NG-PD 15 mm


Let's begin with a relatively basic case (Fig. 2a).

### The initial local injection determines the outcome.

#### Movie 1 : Points in local injection >

Introduction into the submucosal layer is widely regarded as difficult in colorectal ESD. What physicians mean by this is that it can be hard to penetrate the endoscope into the submucosal layer. In this regard, the initial local injection is crucial. If the local injection reliably elevates the lesion and causes the proximal part of the lesion to swell like a bun, then the first incision can open up the optimal layer immediately. When performing the initial


Fig. 1 Ideal local injection



Perform the injection so that the entire lesion is elevated and that the injected solution is pooled immediately below the lesion. The point is to form a shape like a bun.


local injection, focus on creating the "shape" of that injection yourself.

Fig. 2 shows the initial phases of case 1, illustrating how a properly executed local injection allows wide spreading with a single incision and facilitates subsequent introduction into the submucosal layer.




To perform the local injection as shown on the left, inject at a position slightly proximal to the lesion. It is important to hold the mucus gently with the sheath of the injection needle and let the injected liquid flow toward the lesion.


Fig. 2




a: 15 mm LST-NG at rectal Ra.



b: Local injection from the proximal side of lesion to form a bun shape.



c: View of properly performed local injection.



d: Spreading achieved with a single mucosal incision using DualKnife J.

### Proper introduction into the submucosal layer determines the outcome!

#### Movie 2 : Tips for introduction into the submucosal layer and dissection >

The importance of proper local injection is as described above, but proper "introduction into the submucosal layer", i.e., entry of the scope tip beneath the mucosal flap, also leads to a favorable result, after which you simply proceed with dissection (Fig. 3). Merely incising the mucosa does not allow introduction into the submucosal layer, even if you try to insert the hood by force. Introduction into the submucosal layer may be possible after removing submucosal tissue a few times in the optimal manner, but this means that you must hold the proximal mucosa so that a large portion of the submucosal layer which you are going to cut is exposed, as shown in Fig. 4. In this figure, the mucosa is carefully held with the DualKnife J's knob-shaped tip.

A similar effect can be obtained by properly using the 6 o'clock position of the distal hood. Even




Fig. 3: Once the mucosal flap has been securely formed and the endoscope can penetrate into the submucosal layer, subsequent dissection is relatively easy.



Fig. 4: Gently press down the proximal mucosa with the tip of DualKnife J to expose the submucosal layer to be penetrated.

so, it is still not easy to penetrate the layer with a single sweep. Among the techniques that would facilitate introduction into the submucosal layer, one that can be used with DualKnife J is to "lift up with the tip", as shown in Fig. 5. Lift the lesion flap with the sheath of DualKnife J and, after visually recognizing the submucosal layer you are going to cut, adjust the position as if aligning the knob-shaped tip with that line. The key here is to achieve maximum deaeration (introduction into the submucosal layer is impossible if the lumen is completely filled with insufflated air) and to gently lift the layer with the tip of the hood (Fig. 6).



Fig. 5: Insert the sheath of DualKnife J to lift the lesion slightly.




Fig. 6: Repeating incision and dissection a few times enables introduction into the submucosal layer below the mucosal flap. The point is to align the device exit point (forceps channel position) with the submucosal layer to be cut.


### Do not immediately proceed to circumferential incision

The first thing to do after introducing the DualKnife J into the submucosal layer is to excavate the blue submucosal layer in front of you toward the most distal position possible. Rather than immediately proceeding to circumferential incision, dig and advance in the central part immediately beneath the lesion. This is known as the pocket-creation method (PCM)<sup>1</sup>. After widening


the pocket to a certain extent, perform dissection with the goal of opening the pocket to both the left and right (Fig. 7).

What is important here is to advance dissection by making use of the DualKnife J's precision. Your goal should be to remove fibers one by one along the optimal dissection line.


Fig. 7




a: Dig the center of lesion as deeply as possible.



b: After digging the central part for a certain degree, proceed to circumferential incision.



c: Advance dissection with the idea of widening from the center toward the left and right.



d: En-block resection is completed.

## [Case 2] Descending colon IIC 25 mm

When fibrosis is present, introduction into the submucosal layer becomes even more difficult. This case has a higher degree of difficulty than Case 1. It is a IIC lesion near an anastomotic site

of the descending colon after surgery, accompanied by severe cramping and depression (Fig. 8).

### For quick introduction into the submucosal layer – Water-pressure method (WPM)


#### Movie 3 : Water-pressure method >

In this case, introduction into the submucosal layer appears more difficult than usual because the mucosa is thickened due to inflammation. We therefore applied the WPM<sup>2</sup>. This technique was originally proposed for duodenal ESD but can be a powerful tool when introduction into the submucosal layer is difficult. As shown in Fig. 9, it effectively visualizes the tissue to be cut by properly using the tapered hood. Because immersion

expands the visual field, the knob-shaped tip of DualKnife J can be precisely applied to the point to be cut, enabling truly precise ESD.


Once the DualKnife J has been introduced into the submucosal layer, you can use the PCM to excavate the central part and expose the fibrotic region (Fig. 9). After that, you can spread the pocket in the usual way and complete en-bloc resection.

Fig. 8




25 mm descending colon IIC lesion accompanied by severe cramping and depression


Fig. 9




a: The extent of the lesion is hard to identify so peripheral marking is performed.



b: Ideal local injection to form a bun-like shape on the proximal side of the lesion.




c: Quick introduction into the submucosal layer using the WPM with precise application of the DualKnife J tip to the target tissue.



d: Clear exposure of the fibrotic region by using the WPM.



e: Breaching the fibrotic region.



f: En-block resection is complete.