










ORBEYE Room Set-up Guide



SONY 4K 3D
55" Monitor

SONY 4K 3D
55" Monitor

Surgeon

Sterile Table

Assistant


Scrub Nurse


Instruments

Anesthesia

ORBEYE
Video Microscope

CRANIAL: FRONTAL TEMPORAL





CRANIAL: LATERAL SUBOCCIPITAL

SONY 4K 3D
55" Monitor

Assistant

ORBEYE
Video Microscope

Anesthesia


SONY 4K 3D
55" Monitor

Surgeon

Scrub Nurse

Instruments

CRANIAL: MIDLINE SUBOCCIPITAL



SONY 4K 3D
55" Monitor

SONY 4K 3D
55" Monitor

Anesthesia

Assistant

Scrub Nurse

Surgeon


Navigation

ORBEYE
Video Microscope

Sterile Table

Instruments

CRANIOTOMY SETUP



SONY 4K 3D
55" Monitor

Anesthesia

Assistant

ORBEYE
Video Microscope

Instruments

Sterile Table


C-Arm
X-Ray

SONY 4K 3D
55" Monitor


Surgeon

Scrub Nurse

LATERAL MIS




SONY 4K 3D
55" Monitor




Sterile Table



ORBEYE
Video Microscope




Surgeon



SONY 4K 3D
55" Monitor



Assistant



Scrub Nurse



Instruments

PLASTIC FREE FLAP SETUP

SONY 4K 3D
55" Monitor

Scrub Nurse

Assistant

Anesthesia

Instruments

Sterile Table


Surgeon


C-Arm
X-Ray

SONY 4K 3D
55" Monitor


ORBEYE
Video Microscope

SPINE: ANTERIOR CERVICAL





SPINE: LUMBAR MIS & OPEN



UROLOGY MICRO SETUP